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## DX-3100 Series Vibration Sample Magnetometer

### I. Introduction.


DX-3100 series vibration sample magnetometer, is the first set of magnetic moment measurement noise into  $10^{-8}$  emu ( $4 \times 10^{-8}$  emu noise floor at 10 s/point) commercial VSM. We draw lessons from the advantages of several most advanced vibration sample magnetometers in the world, and refer to the design concept of vibration sample magnetometers on superconducting magnets, and design a new generation of products with first-class performance, balanced function and rich selection on the platform of electromagnets.

### II. Functions and characteristics.

1. Hysteresis loop (M-H) can be tested and parameters such as coercivity, residual magnetization and saturation magnetic moment can be obtained.
2. A magnetic moment temperature curve (M-T) can be tested, a superconducting transition temperature can be obtained with a low temperature selection, and a Curie temperature point can be obtained with a high temperature selection.
3. The system is equipped with  $360^\circ$  vibration heads for continuous electric rotation, and the samples are in the center position of the collecting coil and magnetic field during the continuous rotation. Compared with the VSM in the form of superconducting, it is more convenient to study anisotropy.


4. The special design of electromagnet magnetic field gap adjustment mode to ensure that the sample rod in the left and right direction no longer need to adjust the center; extremely precise mechanical device to ensure that the sample rod is in the front and back center of the collection coil each time; the new concept of electric lifting system, can automatically complete the upper and lower direction of the sample center position adjustment; the above three aspects of the design can greatly improve the accuracy and repeatability of the test.


5. The speed of the test has been revolutionized by collecting magnetic moments, magnetic fields and temperature values at 10ms intervals in a row to ensure that users capture any details of the rapid changes in the test process; the rate of change in the electromagnet's magnetic field is as high as 5000 Oe/s, a rate that superconducting magnets can not have; and a hysteresis loop can be obtained at 14 s at this high rate of change in the magnetic field and rate of data acquisition.

6. High precision laser displacement sensor is introduced, which can distinguish the amplitude change of 1 micron, and form closed loop control with the special customized sample vibration drive system, which can achieve 10mm amplitude long time (24 hours) 5/10000 stability.

7. The option replacement and alignment device is designed to ensure that the sample is perpendicular to the ground and does not "hit the wall" when all the options are used, and the vector magnetic moment test is also supported when the high and low temperature option is used.
8. Integrated magnetic moment signal acquisition coil and probe of Gauss meter, less space, more convenient installation, more reliable cable.
9. With a patented design, one-handed sample replacement ease of operation, and software one-click test of the human-computer interface.


### III. Software Interface.



### IV. Specifications and Parameters.

Magnetic Moment Measurement						
Noise	0.1s/pt(RMS)	0.1s/pt(P-P)	1s/pt(RMS)	1s/pt(P-P)	10s/pt(RMS)	10s/pt(P-P)
Room temperature (@3.5mm)	150 nemu	800nemu	80 nemu	380nemu	40 nemu	80nemu
High and low temperature (@8mm)	TBD	TBD	TBD	TBD	TBD	TBD
High or low temperature (@8mm)	0.5 $\mu$ emu	2.5 $\mu$ emu	TBD	TBD	0.1 $\mu$ emu	0.2 $\mu$ emu
Dynamic range	40*10 <sup>-9</sup> emu to 10 <sup>3</sup> emu					
Stability	±0.05% of full scale/day (Fixed coil spacing, data measured 24 hours after constant temperature and magnetic field relative to the maximum range and requiring 2 hours of thermal equilibrium.)					

<b>Repetitiveness</b>	<p style="text-align: center;"><math>\pm 0.5\%</math></p> <p>(The same coil spacing, the same range of repeated loading and removal of samples 10 times measured data.)</p>
<b>Accuracy</b>	<p style="text-align: center;"><math>\pm 1\%</math> (Contrast NIST 772a standard sample)</p>
<b>Sample mass</b>	<p style="text-align: center;">0 to 10g</p> <p>(A larger mass of the sample leads to a decrease in amplitude, which leads to a decrease in SNR.)</p>
<b>Other parameters</b>	<p>System sampling frequency 0.01~10s (continuous mode 10ms);  Time constant of Lock-in amplifier fixed at 100ms, low pass steep drop 18dB/oct;  Vibration amplitude 5~10mm (default 8mm), vibration frequency 55~85 Hz (default 57Hz).</p>

### Magnetic Field Strength

	<b>Magnetic Field Strength</b>			
<b>Sample temperature</b>	<b>Pole cap gap</b>	<b>Sample space</b>	<b>3105(caps dia. 50mm)</b>	<b>3107(caps dia. 50mm)</b>
<b>Room temperature</b>	7.5mm	3.5mm	TBD	3.017T
	12mm	8mm	TBD	2.784T
<b>High and low</b>	20mm	8mm	TBD	2.361T

<b>temperature</b>				
<b>High or low temperature</b>	28mm	8mm	TBD	2.023T

<b>Magnetic Field Performance</b>		<b>Sample Lifting Rotation</b>	
<b>Accuracy</b>	1%	<b>Lifting accuracy</b>	0.1mm
<b>Resolution</b>	1 mOe	<b>Lifting range</b>	0 to 300mm
<b>Stability</b>	100 mOe	<b>Rotation accuracy</b>	0.1°
<b>Rate of change</b>	5000 Oe/s	<b>Rotation range</b>	0 to 730°continuous

<b>Temperature options</b>	<b>Low temperature</b>	<b>High and low temperature</b>	<b>High temperature</b>
<b>Temperature range</b>	80K to 370K	TBD	350K to 1050K
<b>Temperature resolution</b>	0.001K	TBD	0.1K
<b>Stability</b>	±0.1K	TBD	±1K
<b>Time and rate of temperature change</b>	15min min, 6K/min max	TBD	35min min, 20K/min max
<b>Sample environment</b>	Flow nitrogen	TBD	Static argon
<b>Heat-insulating layer</b>	Vacuum	TBD	Vacuum, with multi-layer insulation
<b>Outer diameter and sample space</b>	φ19mm, φ8mm	TBD	φ19mm, φ8mm